

# H<sub>2</sub>O EXPRESS TEST KIT

FAST, EASY, COST-EFFECTIVE SYSTEM RECOMMENDATION

The H<sub>2</sub>O Express Test Kit is designed to help you choose the best Pentair® Everpure® solution for your equipment based on the condition of your water!

Part Number: EV650101 (1PK)  
 Part Number: EV650106 (6 PK)  
 Part Number: EV650206 (6 PK WITH DISPLAY)



Testing is easy! Complete the instruction form included in the kit. Next, collect a water sample according to the instructions and send it to the address on the form. The sample will be tested, and a system recommendation will be sent to the email address(es) provided on the form. Everpure H<sub>2</sub>O Express Test Kits are available through an authorized Pentair Dealer.

## TEST KIT FEATURES:

- X Evaluation of water quality to determine your treatment needs
- X Mail-in lab test for economical analysis and recommendation

## TEST KIT CONTAINS:

- X H<sub>2</sub>O Express Test Kit Form with instructions
- X One (1) 250 mL water sample bottle

Prepaid shipping for Channel partners

## TEST RESULTS INCLUDE:

- |                  |                 |               |             |              |
|------------------|-----------------|---------------|-------------|--------------|
| X pH             | X LSI           | X TDS         | X Chlorides | X Alkalinity |
| X Bicarbonate    | X Calcium       | X Chloramines | X Hardness  | X Sulfates   |
| X Total Chlorine | X Free Chlorine | X Turbidity   | X Silica    |              |

Send your sample to Pentair's test lab.

# HxO Express Test Form

The HxO Express Test Kit contains:

- Customer information and test results form
- One (1) 250mL water sample bottle

## INSTRUCTIONS:

Step 1 - Fill out the section: Test Location, Customer and Dealer Information below.

Step 2 - Fill out the section: Foodservice Equipment below, describing what equipment type the water is being tested.

Step 3 - Place the bottle in the provided sealable bag. Include the customer information, test results form and water sample bottle in the box.

Step 4 - Place the bottle in the provided sealable bag. Include the customer information, test results form and water sample bottle in the box.

Step 5 - Please scan the UPS QR code above for free shipping. If UPS is not available locally, please use the below address: We strongly recommend shipping by a trackable method.

TEST LOCATION AND DEALER INFORMATION			
Test Location Name		City	State Zip Code
Customer Contact Name		Customer Email Address	Customer Phone Number
Dealer Contact Name		Dealer Email Address	Dealer Phone Number

WATER TESTING RESULTS			
Total Dissolved Solids (TDS)		ppm	Silica ppm
Hardness		gpg	pH units
Chloride		ppm	Total - Chlorine ppm
Alkalinity		ppm	Free - Chlorine ppm
Sulfate		ppm	Chloramines (Calculated) ppm
Bicarbonate		ppm	LSI
Calcium		ppm	Turbidity NTU
Pentair Tester Name	Date Received	Lab Test Number	

The purpose of this recommendation is to help Pentair users select the best water treatment for their application based on the water quality parameters presented here and proven technologies to help protect equipment vulnerable to scale, corrosion, and other water quality issues. Pentair is not responsible for equipment conditions at any time. Consult with your Pentair Dealer or contact the factory for further assistance in selecting and sizing the right Pentair solution for your application.



mipl 9EeAiCe... AeUe | (v@3UeA IvA\$e \*4 rImoo | V@eEe^ OeEvEeA | tll>upn>mmqo aVO ?@~a« ñ rol>ol

All indicated Pentair trademarks and logos are property of Pentair. Third party registered and unregistered trademarks and logos are the property of their respective owners.  
©2023 Pentair. All rights reserved.